

IDUC cart

can relabel as 'Urine catheter' if confusion

Stock item

Storeroom code

Comments

Drawer 1



Standard urine catheters: 12, 14, 16

These are slightly too long for the drawers, but work well if the coloured balloon connectors are poking out slightly – this means you can pull out the right size catheter without having to open the drawer which can catch slightly due to the protruding packets. The label on this drawer could also be colour-coded to help with stocking and finding the right size, as it done in the storeroom

There are no special catheters stocked in the trolley – these can be sourced from the storeroom if required e.g. coude tip, three-way, temperature probe etc. These standard catheters are more commonly used and they all contain the correct 10 mL water syringe.

Drawer 2



Standard urine catheters: 18, 20, 22


Notes as above

Drawer 3



| | | |
|--|--|---|
| Catheter tip 60 mL syringes | | |
| Dressing packs | | <i>We have decided to stop using the urine catheter packs as they are oversized and add little beyond standard dressing pack</i> |
| 10 mL water for injection vials | | <i>Do not stock saline as can damage balloon</i> |
| Optilube syringes | | |

Drawer 4



| | | |
|----------------------|--|--|
| Leg bags | | |
| Standard bags | | |
| Urine sample bottles | | |
| 20 mL syringes | | |

Drawer 5



Volumetric hourly urometer reservoir